

NitriSkin XP Nitrile Exam Gloves

12" Powder-Free Non-Sterile Nitrile Examination Gloves

NitriSkin XP Nitrile Exam Gloves is a thicker quality nitrile glove that resists punctures and tears to provide you outstanding barrier protection that is full textured for a more confident grip. When your daily job exposes you to blood-borne pathogens, bacteria and fluids, NitriSkin XP is a glove you can trust.

Characteristics and Overview:

- Full textured for better grip on all items, wet or dry.
- Barrier tested against fentanyl exposure.
- Offer excellent sensitivity for touch perception and dependable protection
- Meet standards of the National Fire Protection Association (NFPA)
- Meets both ASTM D6978-05 and USP 800 Standards for use of chemotherapy gloves when handling hazardous drugs*



Item No	Description	Size	Packaging
MG5008S	NitriSkin XP Nitrile PF Exam Glove	Small	50 eaches/box 10 boxes/case
MG5008M	NitriSkin XP Nitrile PF Exam Glove	Medium	50 eaches/box 10 boxes/case
MG5008L	NitriSkin XP Nitrile PF Exam Glove	Large	50 eaches/box 10 boxes/case
MG5008XL	NitriSkin XP Nitrile PF Exam Glove	X-Large	50 eaches/box 10 boxes/case
MG5008XXL	NitriSkin XP Nitrile PF Exam Glove	2X-Large	45 eaches/box 10 boxes/case
MG5008XXXL	NitriSkin XP Nitrile PF Exam Glove	3X-Large	45 eaches/box 10 boxes/case

Glove Specifications:

Brand	NitriSkin XP
Material	Nitrile
Color	Violet Blue
Sterility	Non-Sterile
Cuff	Beaded Cuff
Sizes	SM, MD, LG, XL, 2XL & 3XL
Accelerators	Yes
Chemo Tested	Yes
Glove Interior	Powder-Free
Glove Exterior	Full Textured
Length (inches)	12 in
Palm Thickness (Avg)	6.1 mil
Finger Thickness (Avg)	8.1 mil
Tensile Strength (Mpa)	38 / 35 Mpa
Ultimate Elongation (%)	650 / 600
Pin Hole (AQL)	1.5

Quality Standards:

Manufactured in accordance with Quality System ISO 9001.
Exceeds and meets current ASTM Standards:
- D6319, D5151, D1671 & D6978.



*Summary of Chemotherapy Testing Reports - ASTM D6978-05

Drug Name	Average Breakthrough Detection Time
Carboplatin	No breakthrough up to 240 min
Cisplatin	No breakthrough up to 240 min
Cyclophosphamide	No breakthrough up to 240 min
Dacarbazine	No breakthrough up to 240 min
Doxorubicin	No breakthrough up to 240 min
Etoposide	No breakthrough up to 240 min
Fluorouracil	No breakthrough up to 240 min
Ifosfamide	No breakthrough up to 240 min
Methotrexate	No breakthrough up to 240 min
Mitomycin C	No breakthrough up to 240 min
Mitoxantrone	No breakthrough up to 240 min
Paclitaxel	No breakthrough up to 240 min
Vincristine Sulfate	No breakthrough up to 240 min
Carmustine	Not for use with Carmustine
ThioTepa	Not for use with ThioTepa