

# NitraGrip Pro Nitrile Exam Gloves

12-inch Powder-Free 2-Ply Blue/White pH Examination Gloves


These 2-Ply technology outer Blue and inner White elongated latex-free nitrile gloves offers the unique ability to detects tears and rips by showing the high visibility white inner through the breach.

## Characteristics and Overview:

- Designed specially for the Emergency Medical Profession for high risk situation that requires the highest standards in double protection and safety in administering critical and patient care.
- Ideal for EMS, Law Enforcement and Corrections
- Meets and exceeds NIOSH recommendations for Fentanyl.
- Meets both ASTM F1342 and D6978-05 Standards for use of chemotherapy gloves when handling hazardous drugs\*
- Meets and exceeds the standard emergency glove medical of NFPA 1999 Standards on protective clothing, Emergency Medical Operations, Edition 2018.



Item No	Description	Size	Packaging
MG50091	NitraGrip Pro Nitrile Exam Glove	Small	50 eaches/box 10 boxes/case
MG50092	NitraGrip Pro Nitrile Exam Glove	Medium	50 eaches/box 10 boxes/case
MG50093	NitraGrip Pro Nitrile Exam Glove	Large	50 eaches/box 10 boxes/case
MG50094	NitraGrip Pro Nitrile Exam Glove	X-Large	50 eaches/box 10 boxes/case
MG50095	NitraGrip Pro Nitrile Exam Glove	2X-Large	45 eaches/box 10 boxes/case
MG50096	NitraGrip Pro Nitrile Exam Glove	3X-Large	45 eaches/box 10 boxes/case



## Glove Specifications:

Brand	NitraGrip Pro
Material	Nitrile
Color	Blue (outer) / White (inner)
Sterility	Non-Sterile
Cuff	Beaded Cuff
Sizes	SM, MD, LG, XL, 2XL & 3XL
Accelerators	Yes
Chemo Tested	Yes
Glove Interior	Powder-Free, pH Balance
Glove Exterior	Full Textured
Length (inches)	12 in
Palm Thickness (Avg)	5.9 mil
Finger Thickness (Avg)	8.0 mil
Tensile Strength (Mpa)	35 / 30 Mpa
Ultimate Elongation (%)	680 / 650
Pin Hole (AQL)	1.5

## \*Summary of Chemotherapy Testing Reports - ASTM D6978-05

Drug Name	Average Breakthrough Detection Time
Carboplatin	No breakthrough up to 240 min
Cisplatin	No breakthrough up to 240 min
Cyclophosphamide	No breakthrough up to 240 min
Dacarbazine	No breakthrough up to 240 min
Doxorubicin	No breakthrough up to 240 min
Etoposide	No breakthrough up to 240 min
Fluorouracil	No breakthrough up to 240 min
Ifosfamide	No breakthrough up to 240 min
Methotrexate	No breakthrough up to 240 min
Mitomycin C	No breakthrough up to 240 min
Mitoxantrone	No breakthrough up to 240 min
Paclitaxel	No breakthrough up to 240 min
Vincristine Sulfate	No breakthrough up to 240 min
Carmustine	Not for use with Carmustine
ThioTepa	Not for use with ThioTepa

## Quality Standards:

Manufactured in accordance with Quality System ISO 9001.  
Exceeds and meets current ASTM Standards:  
- D6319, D5151, D1671 & D6978.

