

NitraPro Nitrile Exam Gloves

16" Powder-Free Non-Sterile Nitrile Examination Gloves

Excellent for sterilization, decontamination, disinfect, cleaning or handling with chemicals in a hazardous environment. This 16-inch specialty glove is designed to offer exceptional added coverage and protection where there is potential exposure to sharps and chemicals.

Characteristics and Overview:


- Full textured fingertips enhance grip in wet and dry situations.
- Longer 16" extended length
- Tested for puncture resistance and color-coded sizing on box.
- Meets both ASTM F1342 and D6978-05 Standards for use of chemotherapy gloves when handling hazardous drugs*



Item No	Description	Size	Packaging
MG50161	NitriSkin Nitrile PF Exam Glove	Small	50 eaches/box 10 boxes/case
MG50162	NitriSkin Nitrile PF Exam Glove	Medium	50 eaches/box 10 boxes/case
MG50163	NitriSkin Nitrile PF Exam Glove	Large	50 eaches/box 10 boxes/case
MG50164	NitriSkin Nitrile PF Exam Glove	X-Large	50 eaches/box 10 boxes/case

Glove Specifications:

Brand	NitraPro
Material	Nitrile
Color	Blue
Sterility	Non-Sterile
Cuff	Beaded Cuff
Sizes	XS, SM, MD, LG, XL and 2XL
Accelerators	Yes
Chemo Tested	Yes
Glove Interior	Powder-Free
Glove Exterior	Full Textured
Length (inches)	16 in
Palm Thickness (Avg)	7.0 mil
Finger Thickness (Avg)	8.9 mil
Tensile Strength (Mpa)	32 / 30 Mpa
Ultimate Elongation (%)	680 / 620
Pin Hole (AQL)	1.5



*Summary of Chemotherapy Testing Reports - ASTM D6978-05

Drug Name	Average Breakthrough Detection Time
Carboplatin	No breakthrough up to 240 min
Cisplatin	No breakthrough up to 240 min
Cyclophosphamide	No breakthrough up to 240 min
Dacarbazine	No breakthrough up to 240 min
Doxorubicin	No breakthrough up to 240 min
Etoposide	No breakthrough up to 240 min
Fluorouracil	No breakthrough up to 240 min
Ifofamide	No breakthrough up to 240 min
Methotrexate	No breakthrough up to 240 min
Mitomycin C	No breakthrough up to 240 min
Mitoxantrone	No breakthrough up to 240 min
Paclitaxel	No breakthrough up to 240 min
Vincristine Sulfate	No breakthrough up to 240 min
Carmustine	Not for use with Carmustine
ThioTepa	Not for use with ThioTepa

